



The “SUITE” Deal Just Got Better!



Thought Technology is proud to launch two new specialized Suites: **EEG Suite** and **Physiology Suite**. The Suites take full advantage of functions and features of the powerful Biograph Infiniti software platform.

They provide complete control over your entire clinical process, including: assessment, training, reporting and demonstrating outcome. For each protocol, the Suites provide signal quality verification screens, multimedia feedback screens and quick report screens showing in-session trends & learning curves.

The EEG Suite contains 6 scripts and more than 100 screens (over 21 from expert clinicians) for 1 or 2-channel EEG on 1 or 2 monitors. It comes with both standard protocols and a variety of template screens for user-defined protocols:

- Standard Bands: Alpha, Alpha peak frequency, Alpha-Theta, Low-high Beta, Beta-Theta, SMR-Theta, Wide-band inhibit, etc.
- 3 User Bands: 1 reward (up), 2 inhibits (down)
- 6 User Bands: 2 rewards (up), 4 inhibits (down)
- Multiple Bands: 1 reward, 9 inhibits

The EEG Suite also offers many advanced tools for assessment and follow-up:

- Pre/Post-training baseline assessment scripts
- A 4-activity EEG assessment script: eyes-open, eyes-closed, and 2 tasks
- Cross-session trend tools to show your client's progress
- A set of hybrid EEG and peripheral integrated biofeedback screens

EEG Suite features include:

- On-the-fly cut-off frequency adjustments
- Impedance detection
- Auto-threshold instruments with user-specified target percentage of time
- Single and dual-monitor display screens
- Real-time epoch statistics for instant trend graphing
- Out-of-this world alpha-theta relaxation script

